

COOLSTORE

New Technology



COOLSTORE 600 / COOLSTORE 1300

Froilabo Firlabo

RELIABLE
STABLE
EFFECTIVE
UNIFORM IN T°

COOLED INCUBATOR

Temperature from 0 up to 60°C



	COOLSTORE 600/0	COOLSTORE 1300/0
Temperature range	from 0°C up to +60°C	
Regulator	PID double display	
Alarms	High and low temperature	
Security	Independent 3.1 overtemperature safety thermostat	
Capacity (l)	600	1350
External dimensions W x D x H (mm)	700 x 830 x 2130	1430 x 830 x 2130
Internal Dimensions W x D x H (mm)	530 x 660 x 1550	1270 x 660 x 1550
Shelves W x D 530 x 655 mm	4	2 x 4
35mm cable port	Right side	
Internal end external cabinet	Stainless steel	
Key lock	Yes	
Castors	4 (2 with breaks)	5 (3 with breaks)
Fluide frigorigène sans CFC / HCFC	R 134a	
Weight (Kg)	170	210
Power supply	1650 W	2400 W
Electric supply	230 V - 50 Hz mono + T	



8 Bd Monge - 69330 Meyzieu Tél: 33 4 78 04 75 75 - Fax: 33 4 78 04 75 76.
Mail : export@froilabo.com

www.froilabo.com

RELIABLE
STABLE
EFFECTIVE
UNIFORM IN T°

COOLED INCUBATOR

Temperature from -20 up to 60°C



COOLSTORE 600/-20

COOLSTORE 1300/-20

Temperature range	from -20°C up to +60°C	
Regulator	PID double display	
Alarms	High and low temperature	
Security	Independent 3.1 overtemperature safety thermostat	
Capacity (l)	600	1350
External dimensions W x D x H (mm)	700 x 830 x 2130	1430 x 830 x 2130
Internal Dimensions W x D x H (mm)	530 x 660 x 1550	1270 x 660 x 1550
Shelves W x D 530 x 655 mm	4	2 x 4
35mm cable port	Right side	
Internal end external cabinet	Stainless steel	
Key lock	Yes	
Castors	4 (2 with breaks)	5 (3 with breaks)
Fluide frigorigène sans CFC / HCFC	R 404a	
Weight (Kg)	170	210
Power supply	1650 W	2400 W
Electric supply	230 V - 50 Hz mono + T	



8 Bd Monge - 69330 Meyzieu Tél: 33 4 78 04 75 75 - Fax: 33 4 78 04 75 76.
Mail : export@froilabo.com

www.froilabo.com