

ARCTIKO

The art of simplicity®

CRYO FREEZER -164/-120°C

ULTF SERIES -164°C

- Digital display
- Key operated
- Password protected operation
- High Temperature alarm
- Low temperature alarm
- Lid open alarm
- Filter alarm
- Contact for remote alarm
- Power failure alarm
- Castors

ULTF1400



SPECIFICATIONS		ULTF 1400
Item Code (DAI no.)		DAI0193
Dimensions exterior WxDxH(mm)		880x1590x1040
Net capacity (L)		128 L
Temperature range (°C)		-120/-164°C
Power supply (V)		380 (3 Phases)
Frequency (Hz)		50

ARCTIKO

The art of simplicity®

Lammefjordsvej 5
DK-6715 Esbjerg N
Denmark

+45 70 20 03 28
+45 70 20 03 29
info@arctiko.com
www.arctiko.com