

# ARCTIKO

The art of simplicity®

## CRYO FREEZER -152/-110°C

### ULTF SERIES -152°C

- Digital display
- Key operated
- Password protected operation
- High Temperature alarm
- Low temperature alarm
- Lid open alarm
- Filter alarm
- Contact for remote alarm
- Power failure alarm
- Castors

### ULTF1500



SPECIFICATIONS		ULTF 1500
Item Code (DAI no.)		DAI0194
Dimensions exterior WxDxH(mm)		890x2005x1009
Net capacity (L)		258 L
Temperature range (°C)		-152/-110°C
Power supply (V)		380 (3 Phases)
Frequency (Hz)		50

## ARCTIKO

The art of simplicity®

Lammefjordsvej 5  
DK-6715 Esbjerg N  
Denmark

+45 70 20 03 28  
+45 70 20 03 29  
info@arctiko.com  
www.arctiko.com